Selective alcohol oxidations with artificial oxidases: alcohol dehydrogenase – NADPH oxidase fusions

Friso S. Aalbers and Marco W. Fraaije
Molecular Enzymology Group
University of Groningen, The Netherlands








ABSTRACT

Selective alcohol oxidations are useful for production of ketones/aldehydes, and for kinetic resolutions of racemic alcohols.[1] Alcohol oxidases (AOXs) and Alcohol Dehydrogenases (ADHs) can catalyze oxidations. Not many AOXs are available, compared to ADHs. ADHs rely on NAD(P)+ for oxidations, and oxidations are thermodynamically less favored – efficient recycling is needed.[1]

Fusion enzymes of three ADHs with an NADPH oxidase (NOX) were produced and characterized (ADHs: *Lb*RADH, *Tb*SADH, ADHA). The fusion enzymes were successfully applied for conversions, and for two activity screening methods: cell extract-based, colony-based.



The fusion enzymes showed good expression with 40-120 mg/L culture (pBAD in *E. coli*). The activity & kinetics were similar to the non-fused enzymes. Conversions were greatly improved by including FAD and catalase.




OUTCOME

- 1. High total turnover numbers (TTN) could be attained for both enzyme and cofactor NADP⁺. Kinetic resolution of *rac*-1-phenyl ethanol could be performed using the NOX-ADH fusions, with high enantiomeric excess.
- 2. Cell extract could be used for rapid activity assay using horseradish peroxidase, without addition of NADP+/NADPH.
- 3. Colonies expressing NOX-ADH fusion could be screened for alcohol oxidation activity active colonies turned dark purple [2].


1. Conversions isolated fusion enzymes


Enzyme	Substrate	X (%)	ee	TTN (enzyme)	TTN (cofactor)
NOX-A	cyclohexanol	95 %	n.a.	31 666	475
NOX-L	cyclohexanol	99 %	n.a.	33 000	495
T-NOX	cyclohexanol	69 %	n.a.	23 000	345
NOX-A	rac-1-phenylethanol	94 %	99% (<i>R</i>)	31 333	470
*NOX-A	rac-1-phenylethanol	56 %	99% (R)	28 000	140
NOX-L	rac-1-phenylethanol	50 %	99% (<i>S</i>)	16 666	250

2. Conversions with cell extracts



3. Facile screening assay with fusion enzymes





CONCLUSIONS

- Expression-level, activity and enantioselectivity were retained.
- The ADH-NOX fusions could be applied as 'artificial oxidases': catalyzing continous alcohol oxidations using only dioxygen, reaching high TTNs.
- With the NOX fusion partner, cell extracts can be used to rapidly measure alcohol oxidation activity through a HRP-coupled assay.
- Colonies expressing ADH-NOX fusion can be screened for specific alcohol oxidation activity – useful tool for screening of large mutant libraries.
- [1] G. Rehn, A. T. Pedersen, J. M. Woodley, *J. Mol. Catal. B Enzym.* 2016, 134, 331–339.
- [2] M. Alexeeva, A. Enright, M. J. Dawson, M. Mahmoudian, N. J. Turner, *Angew. Chemie* 2002, 114, 3309–3312.
- [3] F. S. Aalbers, M. W. Fraaije, submitted, 2018

Conditions

- **1.** Reaction conditions: 5 μ M enzyme fusion, 50 mM substrate, 100 μ M NADP+, 10 μ M FAD, 1 kU catalase, 50 mM Tris/HCl pH 8.5, 0.5 mL volume, 24 °C, 64 hours. *For 24 hours.
- **2.** The reaction included: 50 mM Tris/HCl pH 7.5, HRP (0.8 U), AAP (0.1 mM) and DCHBS (1 mM), cell-free extract containing NOX-A (10 % v/v) and 30 mM cyclohexanol. Controls: no substrate, and no cell-free extract.
- **3.** Method based on [2]. Transformed cells plated onto LB agar plate with 0.2% arabinose and a porous membrane, incubated for 40 hours at 24 °C. Removed membrane with colonies, freeze-thaw step, incubation with assay mix: HRP, 4-chloro-1-naphtol, substrate, NADP+, buffer, 1 % agarose.

"The research for this work has received funding from the European Union (EU) project ROBOX (grant agreement n° 635734) under EU's Horizon 2020 Programme Research and Innovation actions H2020-LEIT BIO-2014-1"

NOTICE: The views and opinions expressed on this poster are only those of the author(s), and do not necessarily reflect those of the European Union Research Agency. The European Union is not liable for any use that may be made of the information contained herein.